
Simple Medical(Shenzhen)Limited Instruction For Use

SIMPLE Knee Brace H4 PCL

Model: 11049& 11050 Size: S, M, L, XL

APPLICATION INSTRUCTIONS

STEP 1. Release the Straps, hold the button 6 and remove the buckle from the buckle seat, and release the strap #4 synchronously. . Release the adjustment belt 7 (see step 1.1) and then release the adjustment seat 8 (the side with the Quick Release Buckle (Fig 1a).

STEP 1.1. Adjust the tightness of the push-back plate. Pull both sides of the rotating knob outward, both sides hear a sound that has unlocked the knob, push the plate back is loosened. Press both sides of the knob inward, both sides hear a sound indicating that the knob has been locked, clockwise of the knob is tightened.

STEP 2. Sit with your legs in the brace, knees slightly bent, center of hinge aligned with center of the knee (Fig 2a).

STEP 3. Tighten the straps. When fastening straps , the product and hinge remain supported to prevent movement. fix the strap #1/#2/#3 in turn. When fastening, the Quick Buckle will be fastened, and the buckle should be fastened in place to ensure the fastening position of the buckle. (Fig 3a) The sequence of tying straps is as follows:

Strap#1: Attach the lower calf position

Strap#2: Attach the upper thigh.

Strap#3: Tie the lower thighs, but not too tight. Excessive tightening may limit circulation or cause the frame to move downward.

STEP 3.1 Adjust the rear push plate. Adjust the seat to the guide rail mouth push back can be fixed. After fixing and locking the knob, rotate the knobs on both sides and adjust the back push plate until the level of the tibial tuberosity (the projecting bone at the upper tibia) is consistent with the level of the healthy leg (adjustment here is for professionals only) (Fig 3b~3c). After adjusting, fix strap #4

STEP 4. Adjust the tightness of each strap for comfort. (Fig 4).

NOTICE:

From the side, the hinge should be slightly behind the midline of the knee. The strap #1, #2 and #4 will affect the tightness of the whole product. In order to ensure good movement, it is normal for the brace to drop during walking. Adjust the position of the straps and bracing if necessary. After the setting is completed, the braces can be taken off by buckles without relaxing and tightening each fixed strap. (Fig 4).

"WEARING TIME"

:Gradually increase the wearing time of the brace and adjust the wearing time of the brace if the body allows.

POSTOPERATIVE WEARING:

The standard of wearing should be determined by the doctor, including the angle of flexion and extension and the time of wearing every day.

The pain caused by osteoarthritis takes time to subside. Wearing a brace regularly for weight-bearing activities (such as walking) will start to feel lighter as it will reduce joint load. This may require wearing it for weeks or more to reduce joint load until your primary care physician allows you to remove the brace.

RANGE OF MOTION ADJUSTMENTS

Start from the standard 0°.Additional flexion and extension angle accessories can be installed and they shall be fitted by the doctor as required. Optional extension angles include 5°, 10°, 15°, 20°, and 25°; Optional flexion angles include 30°, 45°, 60°, 75°, and 90°.

***Please purchase and use the product under the guidance of your doctor or professional.**


IMPORTANT: Both hinges must be set to the same extension and/or flexion settings. For 0° extension limiting, 0° extension stops must be used.

STRAP LENGTH ADJUSTMENT

To determine the required length, open the strap pad, open the end of the binding strap, separate the velcro and the binding strap, cut to the required length, and then attach the strap pad to the end of the binding strap.

TO CHANGE AN EXTENSION STOP:

Unscrew from the front of the hinge cap. Shake out the extension stop. Insert another extension stop into the front of the hinge and replace the screw. Repeat for the opposite hinge. Insert a flexion stop into the back of the hinge and use the extra screw that is attached to the flexion/extension stop tree to secure. Repeat for the opposite hinge.

CARE AND MAINTENANCE

Liners can use neutral detergent and 30 °C water washed by hand, and natural dry. Do not use clothes dryer or other heat sources for drying. Metal parts should not be washed with water. Hinges can be lubricated with dry silicone spray.

WARNINGS:

If you experience pain, swelling, sensory changes, or discomfort with this product, please contact your doctor.

For Single Patient Use Only

NOTICE: The device is a one-element holistic treatment program managed by medical professionals. There is no guarantee that damage can be prevented by using the device.

BASIC PRODUCT ENVIRONMENT

- 1) Ambient temperature: -15 ~ 60 °C
- 2) Relative humidity: ≤ 80%

PRODUCT MATERIALS

- 1)Metal support: aluminum alloy/stainless steel
- 2) Buckle: ABS/ polyester thermoplastic elastomer
- 3) Angle adjustment accessories: nylon + glass fiber
- 4) Strap/pad: nylon/polyurethane/polyester /EVA/ Velcro

CONTRAINDICATION

⚠ Use with caution for unstable fractures;

- ⚠ Use with caution in osteoporosis;
- ⚠ Use with caution for cognitive impairment;
- ⚠ Use with caution for bony ankylosis;


Storage, transportation and expiry date of products

- 1) Term of validity: 5 years.
- 2) Storage of the product: please put it in a dry and clean room before use to ensure indoor ventilation.
- 3) Transportation: the product shall be dustproof, moisture proof, and pollution proof during storage and transportation.

WASTE DISPOSAL

Disposal of wrapping paper (bags): please dispose the Packaging of this product according to the local laws.

SYMBOL AND PICTURE DESCRIPTION

ACCESSORIES LIST

- Instructions ----- 1 copy
- Angle adjusting accessories ----2 PCS (sets)
- Knee pad ---- - 2 PCS (sets)

AFTER-SALES WARRANTY

Warranty: All products from the date of purchase as the original date, normal wear and tear of the product in the use of the process is not considered defective.

The following parts from the date of purchase to the warranty period, if there are any material or process defects, please contact your seller directly. The manufacturer will completely or partially repair or replace the equipment and /or accessories.

Metal support ---- -- 1 year

Simple Medical (Shenzhen) Limited
Address: 519 Wisdom Valley, No.1010 Bulong Road,
Minzhi Street, Longhua District,
Shenzhen, China.518131
Tel : +86-755-86038563
Web: www.simplemedical.cn
E-mail: service@simplegroup.net

Fascinatio Boulevard 522 ,Unit1.7,2909VA Capelle aan den IJssel ,The Netherlands

Tel : +31(0)10-3034500

***Please refer to the Instructions before usage**

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