




SIMPLE
OA KNEE BRACE S1

Low Profile/Superior Quality/Light Weight

+ 3-POINT LEVERAGE SYSTEM

-Simply provide progressive unloading effectively for OA knee patients



+ INDICATIONS

-Mild to severe unicompartmental knee osteoarthritis
-Medial or lateral compartment

+ FEATURES

-Adjustable unloading hinge
-Built-in aerospace-grade aluminum frame
-Adjustable range of motion
-Comfortable non-slip pads
-Unique design
-Latex free & biological evaluation test approved material (ISO10993-5:2009)



Unique Design



Adjustable Unloading Hinge



Adjustable Range of Motion

Comfortable Non-slip Pads



+ ADJUSTMENT PARAMETERS

Varus/valgus angle adjustment range: +/-40°
 Flexion ROM range: 0°, 30°, 45°, 60°, 75°, 90°
 Extension ROM range: 0°, 5°, 15°, 20°, 25°

+ MODELS

Left	11001
Right	11002

+ SIZE

Free Size

+ OA SOLUTIONS



OA Knee Sleeve X1



OA Knee Brace S1



Premium OA Knee Brace S2





Dual Upright
OA Knee Brace M1




SIMPLE OA KNEE BRACE S1

Simple Medical (Shenzhen) Limited

 519 Wisdom Valley, No.1010 Bulong Road, Minzhi Street, Longhua District, Shenzhen, China. 518109

 www.simplemedical.cn

 marketing@simplemedical.cn

 +86-755-86038561