

SIMPLE[®]

Make Healthcare Simple




**SIMPLE
PREMIUM OA
KNEE BRACE S2**

Low Profile/Superior Quality/Light Weight

ISO13485 FDA MDR/CE

+ 3-POINT LEVERAGE SYSTEM

-Simply provide progressive unloading effectively for OA knee patients



+ INDICATIONS

- Mild to moderate unicompartmental osteoarthritis
- Meniscus procedures/repairs
- Medial or lateral compartment

+ FEATURES

- Double hinge adjustment
- One-step closure system
- Built-in aerospace-grade aluminum frame
- Adjustable range of motion
- Non-slip pads
- Unique design
- Latex free & biological evaluation test approved material (ISO10993-5:2009)






Unique Design



One-Step Closure System



Nonslip Pads



Double Hinge Adjustment

+ ADJUSTMENT PARAMETERS

Varus/valgus angle adjustment range: +/-40°
Flexion ROM range: 0°, 30°, 45°, 60°, 75°, 90°
Extension ROM range: 0°, 5°, 15°, 20°, 25°

+ MODELS

Left	11007
Right	11011

+ SIZE

Free Size

+ OA SOLUTIONS



OA Knee Sleeve X1



OA Knee Brace S1



Premium OA Knee Brace S2



Dual Upright
OA Knee Brace M1

SIMPLE PREMIUM OA KNEE BRACE S2

Simple Medical (Shenzhen) Limited

📍 519 Wisdom Valley, No.1010 Bulong Road, Minzhi Street, Longhua District, Shenzhen, China. 518109

🌐 www.simplemedical.cn

✉ marketing@simplemedical.cn

☎ +86-755-86038561