



SIMPLE DUAL UPRIGHT KNEE BRACE H1

Low Profile/Superior Quality

ISO13485 FDA MDR/CE

+ INDICATIONS

- Stabilization of ACL, MCL, PCL,LCL
- Treatment of combined ligament instabilities
- Protection for surgical repairs
- Conservative treatment of ligament damage or deficiency
- Prophylactic bracing for athletic activities


+ FEATURES

- One-step closure system
- Aerospace grade aluminum frame
- Comfortable pads with non-slip
- Adjustable range of motion
0°, 5°, 15°, 20°, 25° extension
0°, 30°, 45°, 60°, 75°, 90° flexion

+ PATENT INFORMATION

- Design Patent NO: 2018305653766
- Utility Model Patent NO: 2018202895414
- PCT NO: PCT/CN2018/079081

+ SIZE



SIZE	THIGH	KNEE	CALF
XS	13~15.5" (33~39cm)	12~13" (33.5~33cm)	10~12" (20.5~30.5cm)
S	15.5~18.5" (39~47cm)	13~14" (33~35.5cm)	12~14" (30.5~35.5cm)
M	18.5~21" (47~53.25cm)	14~15" (35.5~38cm)	14~16" (35.5~40.5cm)
L	21~23.5" (53.25~59.5 cm)	15~17" (38~43 cm)	16~18" (40.5~47cm)
XL	23.5~26.5" (59.5~67.25cm)	17~19" (43~48.25 cm)	18~20" (47~50.75 cm)



Aerospace Grade Aluminum Dual Upright Frame




Nonslip Pads



One-step Closure System



Adjustable Range of Motion



+ ADJUSTMENT PARAMETERS

-Adjustable range of motion
0°, 5°, 15°, 20°, 25° extension
0°, 30°, 45°, 60°, 75°, 90° flexion

+ MODELS

Left	11003
Right	11004

+ SIZE

XS,S,M,L,XL

+ INJURY SOLUTIONS - KNEE BRACE SERIES



Dual Upright Knee Brace H1


Dual Upright Knee Brace H2


Post-Op Knee Brace K1


Post-Op Knee Brace K2


SIMPLE DUAL UPRIGHT KNEE BRACE H1

Simple Medical (Shenzhen) Limited

 519 Wisdom Valley, No.1010 Bulong Road, Minzhi Street, Longhua District, Shenzhen, China. 518109

 www.simplemedical.cn

 joy@simplemedical.cn

 +86-755-86038561