




SIMPLE
OA KNEE SLEEVE X1
Maximum Comfort/Superior Quality

ISO13485 FDA MDR/CE

+ 3-POINT LEVERAGE SYSTEM

-Simply provide progressive unloading effectively for OA knee patients



+ INDICATIONS

- Mild to moderate unicompartmental osteoarthritis
- Meniscus procedures/repairs
- Medial or lateral compartment

+ FEATURES

- Double adjustable hinge even pressure and better comfort
- Comfortable breathable material
- Adjustable range of motion
- Unique varus/valgus angle adjustment, range $\pm 40^\circ$
- Latex free & biological evaluation test approved material (ISO10993-5:2009)





+ ADJUSTMENT PARAMETERS

Varus/valgus angle adjustment range: $\pm 40^\circ$
 Flexion ROM range: $0^\circ, 30^\circ, 45^\circ, 60^\circ, 75^\circ, 90^\circ$
 Extension ROM range: $0^\circ, 5^\circ, 15^\circ, 20^\circ, 25^\circ$

+ MODELS

Left	11005
Right	11006

+ SIZE

S, M, L, XL, XXL

+ OA SOLUTIONS



OA Knee Sleeve X1



OA Knee Brace S1



Premium OA Knee Brace S2



Dual Upright
OA Knee Brace M1



SIMPLE OA KNEE SLEEVE X1

Simple Medical (Shenzhen) Limited

📍 519 Wisdom Valley, No.1010 Bulong Road, Minzhi Street, Longhua District, Shenzhen, China. 518109

🌐 www.simplemedical.cn

✉ marketing@simplemedical.cn

☎ +86-755-86038561